

Filtek™ Supreme Ultra

Ultra Universal Restorative



One simple shading system.
Endless possibilities.

Everyone has a shade, and you've just found a simple way to match it.

As a leader in groundbreaking restorative technologies, 3M ESPE is always in search of ways to improve our products. We've gathered feedback from dentists who use our product every day, to design the next generation of nanocomposite: **Filtek™ Supreme Ultra Universal Restorative.**

If you've never used Filtek Supreme restoratives, you're about to discover what has impressed dentists for years.

If you're already using Filtek™ Supreme Plus Universal Restorative, you'll be happy to know we've made it even better. We've improved the polish retention and enhanced the fluorescence, helping you create natural-looking esthetics.

You'll find Filtek Supreme Ultra restorative more user-friendly—with more single-shade solutions, bold labels that are easier to read, and a wide range of shades and opacities to choose from.

A1B



Simplicity

Exceptionally natural-looking using just one shade.

If you prefer the simplicity of a single-shade technique, you'll be pleased with the unsurpassed esthetics of Filtek™ Supreme Ultra Universal Restorative. Or, if you prefer a multi-shade technique, you'll be equally amazed at how it blends so naturally with your patients' shades.

Each shade of Filtek Supreme Ultra restorative is remarkable by itself or blended with other shades; you'll get the most natural-looking restorations—ever—using your preferred technique.

In addition, we've simplified our shading system. Our syringe and capsule labels are now color-coded by opacity, and the labels have bold lettering to help you quickly find the right shade.

Filtek™ Supreme Ultra

Ultra Universal Restorative

Unsurpassed Esthetics

- Excellent polish and better polish retention than a microfill
- Wide range of shades and opacities
- Improved fluorescence

Simple to Use

- Expanded selection of Body shades
- Color-coded by opacity
- Bold, easy-to-read labels
- Exceptional handling
- Improved handling of Translucent shades

Performance

- Wear-resistance equivalent to Filtek™ Supreme Plus Universal Restorative
- Radiopaque
- Outstanding strength for anterior and posterior use



We've expanded the selection of **Body shades**. Based on dentist requests, we've added darker cervical shades.

Bold new labels are easier to read at a glance, saving you time in finding your selected shade.

We put the **simple** into your restorative technique.



	Body
	A1B
	A2B
	A3B
	A3.5B
	A4B
New	A6B
	B1B
	B2B
	B3B
New	B5B
	C1B
	C2B
	C3B
	D2B
New	D3B
	WB
	XWB

Science

See why 3M ESPE's nanotechnology is so appealing to dentists.


To understand why Filtek Supreme Ultra Universal Restorative delivers results no other composite can, you need to understand HOW our nanotechnology is different, and why it's better.

Nanometer-sized particles are the reason Filtek Supreme Ultra restorative retains its brilliant polish.


The initial gloss of most restoratives can be very impressive. But you want to know your restoration is going to retain its natural gloss for years. Tests prove the polish retention of Filtek Supreme Ultra restorative is better than other composites, and is even better than a microfill.

3M ESPE builds its nanocomposite using a patented process that creates unique clusters of NANOMETER-sized particles. The nanoclusters shear at a rate similar to the wear of


the surrounding resin matrix during abrasion. *This allows the restoration to maintain a smoother surface for long-term polish retention.*



100K magnification of nanocluster



Filtek™ Supreme Ultra nanocluster



Leading competitive composite

SEM photos courtesy of Dr. J. Perdigao, University of Minnesota.

In contrast, most other composites including micro- and nanohybrids incorporate larger particles that contribute significant roughness.




Leading competitive composite

Filtek™ Supreme Ultra Universal Restorative Body shade

Filtek™ Supreme Ultra Universal Restorative Translucent shade

Polish you can see. This side-by-side photo compares the polish retention of a leading restorative to that of Filtek™ Supreme Ultra Universal Restorative after 6,000 cycles of simulated toothbrushing.

Polish retention



Leading competitive composite

Filtek™ Supreme Ultra

This optical analysis compares the surfaces of a leading nano-optimized composite (left) and Filtek™ Supreme Ultra Universal Restorative (Dentin, Enamel, Body shade) after toothbrush abrasion. Source: 3M ESPE internal data

Our “true” nanotechnology provides wear-resistance similar to enamel.

A recent 5-year clinical study* measuring “occlusal contact wear,” showed no statistically significant difference between Filtek™ Supreme Plus Universal Restorative and enamel. Because Filtek Supreme Ultra restorative is derived from the same patented nanotechnology, and in vitro lab data has shown it has a wear-resistance equivalent to Filtek Supreme Plus restorative, *you can be confident using it anywhere in the mouth.*

*Source: Professor Paul Lambrechts and colleagues, Katholieke Universiteit Leuven, Belgium. Data on file.

Versatility

One composite really CAN do it all.

There was a time when you needed a strong composite for posterior fillings, and a highly esthetic composite for anterior work.

Now, ONE composite is strong and esthetic enough to do it all. Imagine how easy it would be to use one composite with complete confidence anywhere in the mouth. ONE composite that is beautiful enough for anterior restorations and strong enough to withstand the rigors of chewing in posterior restorations.

Imagine the possibilities with Filtek Supreme Ultra restorative.

Excellent esthetics and handling led dentists to give Filtek™ Supreme Plus Universal Restorative positive reviews. And the feedback for Filtek Supreme Ultra universal restorative is very positive.

Key improvements include:

- **Remarkable polish:** The restorative is now easier to polish and offers better polish retention than a microfill.
- **Translucent shade handling:** Now our Translucent shades have the same nonsticky, great paste handling as the other opacities—and are also available in capsules.

Find out for yourself how ONE composite outshines, outperforms and outworks, front to back.

Beautifully
versatile
to meet the needs of
your practice.



Color is as much a subjective interpretation as it is a scientific manual on shade determination. Matching does not exist. visual shade. Enamel subjective protocol til others, clinical. a basic knowledge been attained, op. a position to devise individual practice. too often consider the esthetic success. error c might not be noticed. components of color, it is. However, the harmony of the final re in relation to the patient's smile, faci Translucent dental arch as w space, are important. one allows more di perception (e.g., a shade has the sa guide would be a indicate that there discrimination as ch. human observer is dep relative positions of the, the object and the observe

Enamel

Body

Dentin

Opacities: Dentin, Body, Enamel and Translucent designations refer to differences in opacity.



Restoration created using three shades (A3D, A3.5D, A1E) of Filtek™ Supreme Ultra Universal Restorative. Photo courtesy of Dr. Gabriel Krastl, Department of Periodontology, Endodontics and Cariology, University of Basel, Switzerland.



Amalgam filling replaced by Filtek™ Supreme Ultra Universal Restorative in a posterior restoration, creating a strong, esthetic outcome. Photo courtesy of Dr. Grant Chyz.

Ordering Information

Item #	Product Information
6029S	Filtek™ Supreme Ultra Universal Restorative Single Shade Introductory Kit—Capsule Kit includes: 40 – 0.2g capsules: 20 each of Body shades A2 and A3; 1.5g vial Adper™ Easy Bond Self-Etch Adhesive; Technique Guide; Instructions for Use
6028S	Filtek™ Supreme Ultra Universal Restorative Single Shade Introductory Kit—Syringe Kit includes: 2 – 4g syringes: 1 each of Body shades A2 and A3; 1.5g vial Adper Easy Bond Self-Etch Adhesive; Technique Guide; Instructions for Use
6029D Coming soon!	Filtek™ Supreme Ultra Universal Restorative Dual Shade Introductory Kit—Capsule Kit includes: 40 – 0.2g capsules: 10 each of Dentin shades A3, A4, and Enamel shades A2, A3; 1.5g vial Adper Easy Bond Self-Etch Adhesive; 1 Shade Selector Wheel; Technique Guide; Instructions for Use
6028D Coming soon!	Filtek™ Supreme Ultra Universal Restorative Dual Shade Introductory Kit—Syringe Kit includes: 4 – 4g syringes: 1 each of Dentin shades A3, A4, and Enamel shades A2, A3; 1.5g vial Adper Easy Bond Self-Etch Adhesive; 1 Shade Selector Wheel; Technique Guide; Instructions for Use
6029M	Filtek™ Supreme Ultra Universal Restorative Master's Kit—Capsule Kit includes: 129 – 0.2g capsules in plastic reusable case: 3 each of Dentin shades: A1, A2, A3, A4, B3, C4, W 6 each of Body shades: A1, A2, A3, A3.5, A4, B1, B2, C2 3 each of Body shades: A6, B3, B5, C1, C3, D2, D3, W, XW 6 each of Enamel shade: A2 3 each of Enamel shades: A1, A3, B1, B2, D2, W, XW 3 each of Translucent shades: AT, CT 1 Shade Selector Wheel; Technique Guide; Instructions for Use



Syringes and capsule bottle labels are color-coded by opacity for easier shade selection.

Customer Care Center: 1-800-634-2249

www.EveryoneHasAShade.com

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

REFILLS

DENTIN SHADES	SYRINGE	CAPSULE
A1D	6028A1D	6029A1D
A2D	6028A2D	6029A2D
A3D	6028A3D	6029A3D
A4D	6028A4D	6029A4D
B3D	6028B3D	6029B3D
C4D	6028C4D	6029C4D
White	6028WD	6029WD
BODY SHADES	SYRINGE	CAPSULE
A1B	6028A1B	6029A1B
A2B	6028A2B	6029A2B
A3B	6028A3B	6029A3B
A3.5B	6028A3.5B	6029A3.5B
A4B	6028A4B	6029A4B
A6B	6028A6B	6029A6B
B1B	6028B1B	6029B1B
B2B	6028B2B	6029B2B
B3B	6028B3B	6029B3B
B5B	6028B5B	6029B5B
C1B	6028C1B	6029C1B
C2B	6028C2B	6029C2B
C3B	6028C3B	6029C3B
D2B	6028D2B	6029D2B
D3B	6028D3B	6029D3B
White	6028WB	6029WB
Extra White	6028XWB	6029XWB
ENAMEL SHADES	SYRINGE	CAPSULE
A1E	6028A1E	6029A1E
A2E	6028A2E	6029A2E
A3E	6028A3E	6029A3E
B1E	6028B1E	6029B1E
B2E	6028B2E	6029B2E
D2E	6028D2E	6029D2E
White	6028WE	6029WE
Extra White	6028XWE	6029XWE
TRANSLUCENT SHADES	SYRINGE	NOW IN CAPSULES!
Clear	6028CT	6029CT
Blue	6028BT	6029BT
Grey	6028GT	6029GT
Amber	6028AT	6029AT

SYRINGE SHADE REFILLS include 1 – 4g syringe; CAPSULE SHADE REFILLS: Dentin, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use.

Create natural-looking restorations with the easy-to-use 3M™ ESPE™ Shade Selector Wheel. Efficiency is everything in a busy dental practice. Our unique, easy-to-use Shade Selector Wheel simplifies the selection process and serves as a guide for placing one to four shades of composite. To help you meet your patients' needs, we include a Shade Selector Wheel in every Dual-Shade Kit and Master's Kit.



Please recycle. Printed in U.S.A.
© 3M 2010. All rights reserved.
70-2013-0214-1

3M, ESPE, Adper and Filtek are trademarks of 3M or 3M ESPE AG. Used under license in Canada.