PATTERSON® MTA CEMENT - CAPSULES

Endodontic Repair Cement

CIMENT MTA PATTERSON® - CAPSULES

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Instructions for Use

PRODUCT DESCRIPTION

Patterson® MTA Cement is a biocompatible endodontic repair cement in Capsules. Patterson® MTA Cement powder consists of very fine hydrophilic particles of several mineral oxides. After contact with Patterson® MTA Cement liquid, it forms a gel that hardens to an impermeable barrier. Patterson® MTA Cement is delivered in a two-capsule pack. Capsules are easily activated and the content of the capsule is easily ejected out with the Capsule Applicator. Capsule mixing (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/min.

INDICATIONS/INTENDED USE

- Repair of root perforations (during root canal therapy)
- Root end filling (retrograde)
- Pulp capping
- Apexification / Root end filling (orthograde)
- Repair of internal resorption
- Pulpotomy

PERFORMANCE FEATURES

The performance features of the product meet the requirements of the intended use.

CONTRAINDICATIONS

None known

PATIENT TARGET GROUP

Persons who are treated during a dental procedure.

INTENDED USERS

This medical device should only be used by a professionally trained dental practitioner.

SIDE EFFECTS

None known

APPLICATION

1. Activation and Mixing (see Instruction for Capsules)

Activate and mix the Capsule according to the information in the Capsule Cap instruction. Mixing time for the Capsules is 30 seconds.

Attention:

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material. The material must be extruded within 10 seconds after the end of mix.

To prevent dehydration during setting, apply Patterson® MTA Cement intraoral immediately after mixing. Working time (at 23°C) of Patterson® MTA Cement is approx. 2 minutes from start of mix.

2. Applications

2.1. Repair of root perforations

Place rubber dam and clean the root canal system using intracanal instruments and NaOCl solution. Dry the root canal with paper points and isolate the perforation.

Fill the root canal apical of the perforation completely with a suitable root canal filling material.

Mix Patterson® MTA Cement as described under point 1 and extrude it on a glass plate.

Apply Patterson® MTA Cement with suitable instruments into the perforation site and condense it.

Check the position of Patterson® MTA Cement in the root canal by an X-ray. If an adequate barrier has not been created, rinse Patterson® MTA Cement out of the canal and repeat the procedure.

Remove excess moisture with a cotton pellet or a paper point.

Place a cotton pellet in the access to the root canal and apply a temporary filling material. Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done at the earliest 5 minutes after placement of the Patterson® MTA Cement. Patterson® MTA Cement repair material remains as a permanent part of the root canal filling.

2.2. Root End Filling (retrograde)

Create a surgical access to the root-end and resect the root. Prepare an apical cavity to a depth of 3-5 mm.

Isolate the area and dry the root end cavity with paper points. Achieve hemostasis with suitable methods.

Mix Patterson® MTA Cement as described under point 1 and extrude it on a glass plate.

Apply Patterson® MTA Cement with suitable instruments and condense it using a small plugger.

Remove excess cement and clean the surface of the root with a moist piece of gauze.

Confirm placement of the Patterson® MTA Cement repair material with an X-ray. Patterson® MTA Cement remains as a permanent part of the root canal filling.

2.3. Pulp Capping

Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant.

Mix Patterson® MTA Cement as described under point 1 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Patterson® MTA Cement over the exposed pulp and remove excess moisture with a cotton pellet.

At the earliest 5 minutes after application of Patterson® MTA Cement place a small amount of a flowable light cure liner and light cure.

Etch the remaining cavity walls according to the total-etch-technique with an etching gel and apply a suitable bonding agent according to the corresponding instructions.

Place a light cure composite according to the instructions and light cure.

Check pulp vitality and status regularly.

2.4. Apexification / Root End Filling (orthograde)

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal and seal the access opening with a temporary filling material.

After one week remove the calcium hydroxide paste from the root canal.

Mix Patterson® MTA Cement as described under point 1 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Patterson® MTA Cement into the apical region and condense it. Create a 3-5 mm barrier of Patterson® MTA.

Check the position of Patterson® MTA Cement by an X-ray. If an adequate barrier has not been created, rinse Patterson® MTA Cement out of the canal and repeat the procedure.

At the earliest 5 minutes after application of Patterson® MTA Cement place a cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

ADDITIONAL NOTES/WARNINGS

- In the first hour after application handle the placed Patterson® MTA Cement carefully.
- Store Patterson® MTA Cement in the sealed packaging in a dry place prior to use.
- Patterson® MTA Cement can cause discoloration.
- Keep away from children!

COMPOSITION

MTA, bismuth oxide

STORAGE

Store at 10-25°C (50-77°F) in the original seal-pack and protect from moisture. Do not use after expiry date.

Capsules are for single use only.

DISPOSAL

Dispose of the product according to local authority regulations.

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


INSTRUCTION FOR ACTIVATING AND MIXING CAPSULES

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


1.  Capsule avant l'activation.


Poudre
Liquide
Piston


2.  For activation of the Capsule press the plunger on a hard and plane surface to the end into the Capsule.

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

3.  Click once to standardize!
Insert the Capsule into the Capsule Applicator and click once to standardize. Note: The plunger must be at the same level as the bottom of the Capsule.

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

4.  Insert the Capsule into a mixer (or an amalgamator), close lid and mix immediately for 30 seconds (about 4,300 oscillations/min).

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5.  Remove the pin from the nozzle. If not, Capsule can burst.

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6.  Insert the Capsule into the Capsule Applicator. Pull the lever 2 times (2 clicks) to prime the Capsule. Extrude the mixed material on a glass plate and apply directly. Unlock the gun (push button A) and remove the Capsule. Only with the Capsule Applicator the optimal amount of mixed material is guaranteed.

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