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PD 2-05-2014

## Patterson® MTA Cement Endodontic Repair Cement

Patterson® MTA Cement is an endodontic repair cement in a two-capsule pack. Patterson® MTA Cement powder consists of very fine hydrophilic particles of several mineral oxides. After contact with Patterson® MTA Cement liquid it forms a gel that hardens to an impermeable barrier.

Patterson® MTA Cement is delivered in a two-capsule pack. Capsules are easily activated and the content of the capsule is easily ejected out with the application gun. Capsule mixing (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/min.

### 1. Indications

- Repair of root perforations during root canal therapy
- Root-end fillings
- Pulp capping
- Repair of root canals as an apical plug during apexification

### 2. Contraindications

Not known

### 3. Side effects

Not known

### 4. Activation and Mixing (see Instruction for Capsules)

Activate and mix the capsule according to the information in the capsule instruction.

Mixing time for the capsule is 30 seconds.

### Attention:

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting, which may impair or prevent application of the material. The material must be extruded within 10 seconds after the end of mix. To prevent dehydration during setting, apply Patterson® MTA Cement intraorally immediately after mixing.

Working time of Patterson® MTA Cement is approx. 2:00 minutes (at 73°F/23°C).

### 5. Application

#### 5.1. Repair of root perforations

Place rubber dam and clean the root canal system using intra-canal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation.

Fill the apical canal space up to the perforation completely with a suitable root canal filling material.

Mix Patterson® MTA Cement as described under point 4 and extrude it on a glass plate.

Apply Patterson® MTA Cement with suitable instruments into the perforation site and condense it.

Check the position of Patterson® MTA Cement in the root canal by an X-ray. If an adequate barrier has not been created, rinse Patterson® MTA Cement out of the canal and repeat the procedure.

Remove excess moisture with a damp cotton pellet or a paper point.

Place a damp cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the access preparation with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done not before 5 minutes after placement of the Patterson® MTA Cement.

Patterson® MTA Cement repair material remains as a permanent part of the root canal filling.

#### 5.2. Root-End Filling

Create an access to the root-end and resect the root with a surgical bur.

Use an ultrasonic tip to prepare a class I root-end cavity preparation to a depth of 3-5 mm.

Isolate the area and dry the root end cavity with paper points.

Patterson® MTA Cement as described under point 4 and extrude it on a glass plate.

Apply Patterson® MTA Cement with suitable instruments and condense it using a small plugger.

Remove excess cement and clean the surface of the root with a moist piece of gauze.

Confirm placement of the Patterson® MTA Cement repair material with an X-ray.

The Patterson® MTA Cement repair material remains as a permanent part of the root canal filling.

## Directions for Use

### 5.3. Pulp Capping

Place rubber dam and prepare the cavity outline. If caries is present, remove it. Rinse cavity and exposed pulpal areas with a suitable disinfectant.

Mix Patterson® MTA Cement as described under point 4 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Patterson® MTA Cement over the exposed pulp and remove excess moisture with a dry cotton pellet.

Not before 5 minutes after application of Patterson® MTA Cement place a small amount of a flowable light cure liner and light cure.

Etch the remaining cavity walls according to the total-etch-technique with an etching gel and apply a suitable bonding agent according to the corresponding instructions.

Place a light cure composite according to the instructions and light cure.

Pulp vitality and status should be checked by X-ray at regular intervals.

### 5.4. Apexification

Place rubber dam and clean the root canal system using intra-canal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal for one week. Seal the access opening with a temporary filling material.

Mix Patterson® MTA Cement as described under point 4 and extrude it on a glass slab.

With a suitable instrument apply a small amount of Patterson® MTA Cement into the perforation site and condense it. Create a 3 - 5 mm barrier of Patterson® MTA Cement.

Check the position of Patterson® MTA Cement by an X-ray. If an adequate barrier has not been created, rinse Patterson® MTA Cement out of the canal and repeat the procedure.

Remove excess moisture with a damp cotton pellet or a paper point.

Place a damp cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the access preparation with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done not before 5 minutes after placement of the Patterson® MTA Cement.

Patterson® MTA Cement repair material remains as a permanent part of the root canal filling.

### Additional remarks

• Store Patterson® MTA Cement in the sealed packaging at a dry place prior to use.

• Intraoral application of Patterson® MTA Cement must be done immediately after mixing to prevent dehydration during setting.

• Patterson® MTA Cement can cause discoloration. Use Patterson® MTA Cement only in the root canal and/or the pulp chamber.

• In order to obtain a creamy consistency, 1 or 2 drops of sterile water can be added to the mixed Patterson® MTA Cement.

### Storage

Store Patterson® MTA Cement at a dry place at 50-77°F/10-25°C. Do not store below 50°F/10°C!

Do not use after expiration date.

Capsules are for single use only.

### Keep away from children!

For dental use only!

### Caution:

Federal law restricts the sale of this device to or by the order of a dentist.

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







## Instruction for activating and mixing of capsules

Capsules before activation.


For activation of the capsules, press the plunger on a hard and plane surface to the end into the capsules.


Insert the capsules into the application gun and **click once** to standardize.


**Note:** The plunger must be at the same level as the bottom of the capsule.

Insert the capsules into the amalgamator, close lid and mix immediately for 30 seconds (about 4300 oscillations / min.)

Remove the pin from the nozzle. If pin is not removed, capsule can burst.

Insert the capsules into the application gun. Pull the lever 2 times (2 clicks) to prime the capsules. Extrude the mixed material directly into the preparation. Unlock the gun (push button A) and remove the capsules.

**Only** with the application gun is the optimal amount of mixed material is guaranteed.

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Patterson Dental Supply, Inc.