Polyvinylsiloxane Impression Material Cartridge Delivery System

DIRECTIONS FOR USE (Addition Reaction, Type 1)

INTRODUCTION

PATTERSON BRAND Bite Registration addition cured silicone elastomeric impression material meets or exceeds all requirements for Type 1 impression materials specified by the American Dental Association Specification No. 19.

Addition materials are manufactured to perform according to their specifications for work and set time at 73° F (23° C). Higher temperature will decrease the work and set time. Conversely, materials that are lower than 73° F will have extended work/set times. The shelf life of all **PATTERSON BRAND** impression materials is two years. During hot summer months, they can be stored in the refrigerator to prevent degradation - but should be allowed to warm up to room temperature.

In any technique where an impression tray is required the first step is to paint the tissue side of the preselected tray with vinyl adhesive. The optimum bond occurs when the adhesive film is as thin and even as possible. The adhesive should be allowed to dry for 5 minutes.

DIRECTIONS FOR BITE REGISTRATION

PATTERSON BRAND Bite Registration is a unique mousse-like impression material that was originally designed for taking bite relationships between the upper and lower jaws. In the profession, it has found much broader use as a general impression material.

This type is the *regular* two minute work time material. If only bite registrations are to be taken, there is a fast version with a 30-second work time and a one minute mouth removal. This type is too fast as a general impression material.

NOTES:

Vinyl impression materials will not work properly and may never set when sulfur, eugenol or amines are present. For example, latex examination gloves - if placed in contact with the teeth or soft tissue will probably inhibit the set of the vinyl. Certain gingival retraction cords which contain sulfides can cause sulcus areas to remain unset.

All vinyl adhesives are flammable and contain, toxic, volatile materials. DO NOT inhale and DO NOT have any open flames near the adhesive during application to the tray. Always close the adhesive bottle immediately after use.

PATTERSON BRAND impression material is completely non-toxic. However, silicone can stain clothing and when set is chemically resistant. Before setting, spots can be removed with commercial dry cleaning fluids.

WORKING AND SET TIMES FOR PRODUCTS

PATTERSON BRAND Bite Registration Regular Set

Work: 2:00 Work/Set time: 3:00

PATTERSON BRAND Bite Registration Fast Set

Work: 0:30 Work/Set time: 1:30

PATTERSON BRAND Bite Registration Super Fast Set

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