3M ESPE

Ultimate simplicity in temporization.







- 9 preformed sizes
- Custom fit
- Excellent strength
- Natural-looking esthetics
- Light curable
- Low intra-oral heat during fitting process
- Low oxygen-inhibited layer
- No unpleasant odor
- Easy to add-on with flowable composite

The world's first preformed, malleable temporary crown.

3M ESPE has brought many firsts to dental professionals worldwide. Now, we bring you Protemp[™] Crown Temporization Material, the world's first single-unit self-supporting, malleable, light-curable composite crown. This revolutionary, time-saving technique can result in higher productivity for you.

As simple as $1-2-3^*$

No impression or matrix. No dispensing hardware or mix tips.



Trim following gingival contour





Finish and polish

What dentists are saying.

G Making temporaries is now very convenient with no preliminary impression. **This revolutionary new process is simple, quick and easy, producing a temporary that fits well, is kind to the tissue,** does not come off in my experience, and keeps the occlusion in place. A definite winner for 3M ESPE! **L** This is a paradigm shift in temporization. Adjectives alone cannot convey the practical application of a product like this. It's a total success for high-quality temporary restorations. **JJ** — Robert G. Ritter, DMD

*Refer to product instructions for complete step-by-step information.

Compare what you can get in 4 minutes or less.

	Protemp [™] Crown	Bis-Acrylic Materials	Powder-Liquid Systems			
Single unit under 4 minutes	V ,	Highly variable; dependent on matrix/impression & temporization material used				
No matrix/impression needed		4	4			
Custom fit						
Preformed anatomy						
Low oxygen inhibition layer		4	4			
Esthetic results	V .					
Low intra-oral heat						
Cure on demand		-				

Ordering Information

	Item #	Product Information					
	50600	Protemp [™] Crown Temporization Material Introductory Kit Kit Includes: 42 Temporary Crowns; Crown Size Tools (2 bags of 10); Instructions; Technique Guide; Crown Drawer; Crown Drawer Labels					
		Temporary Crown Sizes Molar Upper Large - 6 Molar Upper Small - 6 Molar Lower Large - 6 Molar Lower Small - 6	Bicuspid (Premolar) Bicuspid (Premolar) Bicuspid (Premolar)	Upper Small - 4	Cuspid (Canine) Large - 3 Cuspid (Canine) Small - 3		
3M™ ESPE™ Protemp™ Crown Temporization Material Introductory Kit (50600)	50600TK	K Protemp [™] Crown Temporization Material Trial Kit Includes: 16 Temporary Crowns; Crown Size Tools (1 bag of 10); Instructions; Technique Guide; CD					
	Protemp™	Temporary Crown Sizes Molar Upper Large - 3 Bicuspid (Premolar) Upper Large - 2 Molar Upper Small - 3 Bicuspid (Premolar) Upper Small - 2 Molar Lower Large - 2 Bicuspid (Premolar) Upper Small - 2 Molar Lower Small - 2 Bicuspid (Premolar) Lower - 2 Molar Lower Small - 2 Bicuspid (Premolar) Lower - 2 More Temporization Material — Refills					
	50610	Molar Upper Large Kit includes: 5 Temporary Crov Instructions; Crown Size Tool	vns;	50614	Bicuspid (Premolar) Upper Large Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		
3M™ ESPE™ Protemp™ Crown Temporization Material Trial Kit (50600TK)	50611	Molar Upper Small Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		50615	Bicuspid (Premolar) Upper Small Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		
	50612	Molar Lower Large Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		50616	Bicuspid (Premolar) Lower Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		
	50613	Molar Lower Small Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		50617	Cuspid (Canine) Large Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		
				50618	Cuspid (Canine) Small Refill Kit includes: 5 Temporary Crowns; Instructions; Crown Size Tool		
	Protemp [™] 3 Garant [™] Temporization Material						
	46952	Starter Pack Includes: 1 Protemp 3 Garant 50ml (67g) cartridge, A2; 15 Garant 2 mixing tips, blue; 1 Garant 2 Dispenser, complete 4:1/10:1: 1 Protemp 3 Add-On (2n) dispenser A2: BelyX™ Temp Cement NE					

3M[™] ESPE[™] Protemp[™] 3 Garant[™] Temporization Material Starter Pack (46952)

complete 4:1/10:1; 1 Protemp 3 Add-On (2g) dispenser, A2; RelyX™ Temp Cement NE

Customer Care Center: 1-800-634-2249

www.3MESPE.com/innovation



Dental Products

3M Center Building 275-2SE-03 St. Paul, MN 55144-1000 U.S.A. 1-800-634-2249

3M Canada Post Office Box 5757

London, Ontario N6A 4T1 Canada 1-888-363-3685

Please Recycle

Printed in U.S.A. © 3M 2007. All rights reserved. 3M, ESPE, Garant, Protemp and RelyX are trademarks of 3M or 3M ESPE AG. Used under license in Canada.

70-2009-3866-3