


Have You Experienced Operative Carbide Breakage? **WELDED CARBIDES ARE STRONGER!**




Instruments above are shown larger than actual size.

NEW & IMPROVED NTI 2-Piece Operative Carbides are welded, not soldered, and are stronger at the joint between the cutting head and stainless-steel shank resulting in a smoother, more precise cutting action. See for yourself why NTI 2-Piece Operative Carbides differ from others who use soldering technology! Available in Friction Grip, Friction Grip Short-Shank, Right Angle/Latch, Surgical Length FG, Surgical Length RA, and Handpiece shanks. Friction Grip, Friction Grip Short-Shank, and Right Angle/Latch are available in 5-packs or Big Packs (100 of same) in selected shapes.


What is a "Sharp" version?

The "Sharp" version NTI 2-Piece Operative Carbides contain a transversing blade on the tip of the instrument. The transversing blade is designed for superior end-cutting, sectioning tooth structure, and endodontic access. The "Sharp" version is available on round shapes and in Right Angle/Latch shank only.

Top View of H10S-RA (Sharp)



Top View of H10-RA



	ROUND											
U.S. NO.	1/4	1/2	1	2	2 Sharp	3	3 Sharp	4	4 Sharp	5	5 Sharp	6
SIZE 1/10mm	005	006	008	010	010	012	012	014	014	016	016	018
LENGTH mm	-	-	-	-	-	-	-	-	-	-	-	-
FG (5-pk)	H0.25-FG	H0.5-FG	H1-FG	H2-FG				H4-FG				H6-FG
FG (100-pk)	H0.25-FG100	H0.5-FG100	H1-FG100	H2-FG100				H4-FG100				H6-FG100
FG SS (5-pk)				H2-FGSS				H4-FGSS				H6-FGSS
FG SS (100-pk)								H4-FGSS100				H6-FGSS100
FG SG (5-pk)			H1-FGSG	H2-FGSG				H4-FGSG				H6-FGSG
RA (5-pk)	H0.25-RA	H0.5-RA	H1-RA	H2-RA	H2S-RA	H3-RA	H3S-RA	H4-RA	H4S-RA	H5-RA	H5S-RA	H6-RA
RA (100-pk)				H2-RA100				H4-RA100				H6-RA100
RA SG (5-pk)				H2-RASG				H4-RASG				H6-RASG
HP (5-pk)		H1-006HP		H1-010HP				H1-014HP				H1-018HP

	ROUND						INVERTED CONE				PEAR	
U.S. NO.	7	7 Sharp	8	8 Sharp	10	10 Sharp	33 1/2	34	35	37	330	331
SIZE 1/10mm	021	006	023	023	027	027	006	008	010	014	008	010
LENGTH mm	-	-	-	-	-	-	0.6	0.8	1.0	1.4	1.6	1.7
FG (5-pk)			H8-FG				H33.5-FG	H34-FG	H35-FG	H37-FG	H330-FG	H331-FG
FG (100-pk)			H8-FG100				H33.5-FG100	H34-FG100	H35-FG100		H330-FG100	H331-FG100
FG SS (5-pk)									H35-FGSS		H330-FGSS	
FG SS (100-pk)									H35-FGSS100		H330-FGSS100	
FG SG (5-pk)	H7-RA	H7S-RA	H8-FGSG	H8S-RA	H10-RA	H10S-RA	H33.5-RA	H34-RA	H35-RA	H37-RA		
RA (5-pk)			H8-RA									
RA (100-pk)			H8-RA100					H2-008HP	H2-010HP	H2-014HP	H7-008HP	
RA SG (5-pk)			H8-RASG									H7-012HP
HP (5-pk)			H1-023HP									

	LONG PEAR	CYLINDER		CROSS CUT CYLINDER				ROUND END CYLINDER				
U.S. NO.	245	56	57	58	556	557	558	559	1156	1157	1158	1159
SIZE 1/10mm	008	009	010	012	009	010	012	014	008	010	012	014
LENGTH mm	2.0	4.2	4.2	4.2	4.2	4.2	4.2	4.2	4.2	4.2	4.2	4.6
FG (5-pk)	H245-FG	H56-FG	H57-FG	H58-FG	H556-FG	H557-FG	H558-FG		H1156-FG	H1157-FG	H1158-FG	H1159-FG
FG (100-pk)	H245-FG100	H56-FG100		H58-FG100	H556-FG100	H557-FG100	H558-FG100		H1156-FG100	H1157-FG100	H1158-FG100	H1159-FG100
FG SS (5-pk)	H245-FGSS	H56-FGSS			H556-FGSS	H557-FGSS						
FG SS (100-pk)	H245-FGSS100	H56-FGSS100			H556-FGSS100	H557-FGSS100						
FG SG (5-pk)						H557-FGSG	H558-FGSG	H559-FGSG				H1157-FGSG
RA Latch (5-pk)						H557-RA						
RA SG (5-pk)						H557-RASG	H558-RASG					H1157-RASG
HP (5-pk)					H31-009HP	H31-010HP	H31-012HP					H21R-014HP

	TAPERED FISSURE											
U.S. NO.	1558	168	169	169L	170	171L	699	700	701	701L	702	703
SIZE 1/10mm	012	006	009	009	010	012	009	010	012	012	016	021
LENGTH mm	4.2	4.2	4.2	5.2	4.2	6.0	4.2	4.2	4.2	6.0	4.6	4.9
FG (5-pk)	H1558-FG		H169-FG	H169L-FG			H699-FG	H700-FG	H701-FG		H702-FG	
FG (100-pk)	H1558-FG100		H169-FG100	H169L-FG100					H701-FG100			
FG SS (5-pk)									H701-FGSS			
FG SS (100-pk)												
FG SG (5-pk)	H1558-FGSG							H700-FGSG	H701-FGSG		H702-FGSG	H703-FGSG
RA Latch (5-pk)		H168-RA			H170-RA							
RA SG (5-pk)												
HP (5-pk)	H31R-012HP		H23-009HP		H23-010HP	H23L-012HP		H33-010HP	H33-012HP	H33L-012HP	H33-016HP	H33-021HP

FG: Friction Grip FG SS: Friction Grip Short Shank FG SG: Friction Grip Surgical RA: Right Angle/Latch RA SG: Right Angle/Latch Surgical HP: Handpiece

*FREE GOODS MUST BE OF EQUAL OR LESSER VALUE. FREE GOOD OFFERS CANNOT BE COMBINED WITH ANY OTHER OFFERS. QUALIFYING ORDERS MUST BE ON ONE INVOICE; MULTIPLE INVOICES CANNOT BE COMBINED. To receive free goods, fax or mail a copy of your preferred authorized dealer invoice dated between January 1, 2005 and August 31, 2005, along with your name, address, phone number and email to AXIS Dental (Free Goods Offer), 3008 W. Story Rd., Irving, TX 75038 or FAX to (972) 257-3647. All copies must be received no later than September 15, 2005. Allow 6-8 weeks for delivery. All free good offers are applicable to retail customers only. PLEASE CONTACT YOUR PREFERRED AUTHORIZED DEALER FOR A COMPLETE SELECTION OF ALL AXIS/NTI PRODUCTS.