

2017 CLINICAL PROBLEM SOLVER

Maxcem Elite™ Chroma

Kavo Kerr



www.kerrdental.com

Clinician: Dr. Lori K. Brown



CLINICAL PROBLEM:

It's hard to know when the ideal time is to start removing excess cement after seating a restoration. If you wait too long, clean up can be a nightmare.

SOLUTION:

Maxcem Elite™ Chroma is the first self-adhesive resin cement with color changing chemistry that lets you know when the cement is ready for clean up. The cement dispenses pink but fades in color when it reaches its gel state, indicating the most ideal time to remove excess material. Once it is ready to remove, the excess cement peels off in one motion. This results in less discomfort to the patient and saves time... a game changer!