

## What's Hot and What's Getting Hotter

BY HOWARD S. GLAZER, DDS, FAGD

### MATERIAL

#### **Flexitime® Bite**

Heraeus Kulzer, LLC  
300 Heraeus Way  
South Bend, IN 46614  
800.431.1785

[www.heraeus-kulzer.com](http://www.heraeus-kulzer.com)

There have been many vinyl polysiloxane (VPS) bite registrations to come to market over the past few years. They are extremely useful and necessary and have become faster and better over time. However, this is the first bite registration material I am aware of that offers superior powder-free scannability and optimization



for use with computer-aided design (CAD) technology.

Flexitime® Bite, of course, comes in a "no fuss no muss" automix cartridge that fits the standard 1:1/2:1 mixing gun. Flexitime Bite can be used with Flexitime and Flexitime Xtreme, as well as most other VPS materials. It has a working time of 30 seconds and fully sets in 30 seconds. It registers an impressive durometer of 40, which

means it will not yield and can be trimmed and replaced on and off a working model many times. Flexitime Bite is purple in color, allowing for easy portrayal of the opposing dentition so that

accuracy of the bite is maintained. As with any impression, it should be placed in a disinfection solution, such as 5.25 percent sodium hypochlorite or similar, for 10 minutes before being sent to the laboratory. Accuracy of the occlusion interface is essential for proper function and patient comfort. These features make the Flexitime Bite registration material quick and easy to use.

# Heraeus



## Flexitime® Bite

Powderless scanning for use with CEREC® Technology.

**Flexitime Bite** delivers exceptional accuracy for standard bite registrations and superior powder-free scannability for use with CAD technology.

- Extremely high shore D hardness of 40
- Sets fast in just 30 seconds
- Whipped cream-like consistency ensures exceptional accuracy
- Easy to trim and does not crumble

*"Flexitime Bite can be used without any dusting for a predictable result."*

**Sameer Puri, DDS**  
CERECdoctors.com Co-Founder



## Flexitime®

Flexitime® and Venus® are trademarks of Heraeus Kulzer GmbH. CEREC® is a trademark of Sirona.

Heraeus Kulzer, LLC 300 Heraeus Way, South Bend, IN 46614 Phone: (800) 431-1785 Fax: (877) 271-5211 [www.heraeus-kulzer-us.com](http://www.heraeus-kulzer-us.com) [www.smilebyvenus.com](http://www.smilebyvenus.com)