

Heraeus

# The details you need.



# The comfort your patients want.

Soft edges and high-contrast make AGFA the patient's choice for comfort and the dentist's choice for diagnosis.

- AGFA Dentus® intraoral film is soft, flexible, and features a radius edge to ensure patient comfort
- High contrast and maximum density with excellent grey graduation
- Intraoral films are fully protected from moisture contamination and can be disinfected



**AGFA**

Patient Preferred.

## Intraoral Film

### PATIENT PREFERRED

- Soft, flexible film features a radiused edge to ensure patient comfort
- Exclusively made with Softopac®, smooth edges and rounded corners eliminate irritation

### OPTIMAL IMAGE SHARPNESS

- High contrast and high maximum density support the diagnosis due to excellent grey graduation
- Flexible but firm films help to reduce distortion

### CONVENIENT

- Double sealed packets are fully protected from moisture contamination
- Can to be immersed in disinfection liquid without worry
- Full offering to meet all your needs—adult PA (size 2), pedo (size 0), intraoral E/F and D Speed, extraoral, duplicating, and accessories

## Extraoral Film

- Very high definition due to disc grain technique
- Very low Base-Plus-Fog Density
- Particular contrasting in various density ranges

### AGFA Dentus M2 Comfort



Speed F/E 0.10 sec.



Speed D 0.20 sec.

### Competitor Film



Speed F/E 0.10 sec.



Speed D 0.20 sec.

Prepared under identical condition to yield optimal density.  
X-ray equipment: Shortcon, 70kV – 7mA, head phantom – mandibular, filter 10 mm.  
Automatic processor: 28°C, 6 minutes.



#### M2 COMFORT — INTRAORAL E SPEED

*Less Radiation! Shorter exposure time with E/F-Speed film*  
**SIZE 0**

**65031124** M2-54 1 Film Per Packet (100 Count)

**65031256** M2-53 2 Films Per Packet (100 Count)

#### **SIZE 2**

**65414615** M2-58 1 Film Per Packet (150 Count)  
Case Lot 24

**65414623** M2-57 2 Films Per Packet (150 Count)

#### **SIZE 4 (OCCLUSAL)**

**65870955** M2-50 1 Film Per Packet (25 Count)

#### M2 COMFORT — INTRAORAL D SPEED

*Traditional exposure time. No change to current technique or settings*

**66009814** M2-D01 AGFA Dentus M2 Comfort D  
(2 x 3) 100 x 1 Sz-0

**66009815** M2-D02 AGFA Dentus M2 Comfort D  
(2 x 3) 100 x 2 Sz-0

**66009817** M2-D1 AGFA Dentus M2 Comfort D  
(3 x 4) 150 x 1 Sz-2 (Case Lot 24)

**66009818** M2-D2 AGFA Dentus M2 Comfort D  
(3 x 4) 150 x 2 Sz-2 (Case Lot 24)

#### PANORAMIC

**65780689** RP6 (5 x 12 in)  
100 Films, Blue Sensitive

**65780697** RP6 (6 x 12 in) (15 x 30 cm)  
100 Films, Blue Sensitive

**66000311** Ortholux (5 x 12 in)  
100 Films, Green Sensitive

**66000312** Ortholux (6 x 12 in) (15 x 30 cm)  
100 Films, Green Sensitive

#### AGFA DUPLICATING FILM

**65038641** AGFA Dentus Duplicating Film  
(5 x 12 in, 100 Films)

**65038668** AGFA Dentus Duplicating Film  
(6 x 12 in, 100 Films)

#### CEPHALOMETRIC

**66000314** Ortholux (8 x 10 in) 100 Films,  
Green Sensitive

#### AGFA INTENSIFYING SCREENS (DENTUS GREEN)

**65769847** Ortho Regular (15 x 30 cm) (6 x 12 in)

#### ACCESSORIES

**66002011** AGFA Dentus Ortholux Cassette (15 x 30 cm)

**66002014** AGFA Dentus Ortholux Cassette (8 x 10 in)

**To learn more about AGFA call 1-877-258-5937**