

Filtek™ Bulk Fill

Posterior Restorative



Fast and easy

One increment up to 5 mm.

No expensive dispensing devices
or time-consuming layers.

Posterior restorations made simple.




Filtek™ Bulk Fill Posterior Restorative

There's no question that bulk fill materials save time, however, true simplification involves eliminating complicated expensive dispensing devices, additional layers and multiple steps, without compromising wear resistance or stress relief.

Filtek™ Bulk Fill Posterior Restorative is your one-step bulk placement solution designed to improve productivity by allowing fast and easy posterior restorations.

- ✓ One-step placement, no additional capping layer
- ✓ Fast and easy procedure
- ✓ Excellent adaptation without additional expensive dispensing devices
- ✓ Stress relief to enable up to 5 mm depth of cure
- ✓ Better in vitro wear resistance than market-leading bulk fill materials
- ✓ Excellent handling and sculptability
- ✓ TRUE nanofiller technology and two innovative methacrylate monomers act to lower polymerization stress without compromising wear

It's as simple as “one and done.”



One Bulk Fill Placement
up to 5 mm




10 sec. light cure x 3
(occlusal, buccal, lingual)

Dentists were asked to place and cure restorations and rate the materials based on preference.

In a simulated operatory setting, we put Filtek™ Bulk Fill Posterior Restorative head-to-head against several restorative materials that require different placement techniques. The materials included an incrementally placed composite and a flowable bulk fill composite.

Filtek™ Bulk Fill Posterior Restorative stands out over the competition.

24%
faster
placement time
vs incrementally
placed composite





2X
number of
dentists

prefer Filtek™ Bulk Fill
Posterior Restorative

vs a bulk fill flowable based
on placement technique

TRUE nanofiller technology and innovative chemistry

5K magnification
of nanocluster



50K magnification
of nanocluster

Filtek Bulk Fill Posterior Restorative uses the same TRUE nanofiller technology that's been time-proven by other leading Filtek™ Restorative products. Not only does this technology contribute to high strength, excellent handling and better wear resistance, but also to restorations that are fast and easy to polish.

In addition to this TRUE nanofiller technology, Filtek Bulk Fill Posterior Restorative also contains two innovative **methacrylate monomers** that, in combination, act to lower polymerization stress, without compromising wear.

The combination of these two technologies results in an easy-to-use, one-step placement, restorative material designed to simplify your posterior restorations.

Also from the brand you trust ...

Filtek™ Bulk Fill

Flowable Restorative

An excellent low-stress liner/base solution.

If you're looking for a great, all-around, low-stress bulk fill flowable, this is it. Filtek™ Bulk Fill Flowable Restorative is up to 50% stronger than the leading bulk fill flowable¹—and offers over 3X the wear resistance.¹ Excellent for use as²:

- a liner for Class I and Class II restorations
- a fast bulk fill base up to 4 mm for Class I and Class II restorations
- as a wear-resistant solution for small, non stress-bearing occlusal surfaces, Class III and Class V

1. 3M ESPE in vitro data on file.

2. Refer to Instructions for Use for complete list of indications.



Ordering Information

Filtek™ Bulk Fill Posterior Restorative

Item #	Product Information
038-7035	Filtek™ Bulk Fill Posterior Restorative Capsule Procedure Kit Kit includes: 40 – 0.2g capsules: 20 each of A2 and A3; 5.0ml vial Scotchbond™ Universal Adhesive; Instructions for Use; 10 – Sof-Lex™ Spiral Finishing Wheels; 10 – Sof-Lex™ Spiral Polishing Wheels; 1 – RA Mandrel; Technique Guide; Instructions for Use
038-7043	Filtek™ Bulk Fill Posterior Restorative Syringe Procedure Kit Kit includes: 2 – 4g syringes: 1 each of A2 and A3; 5.0ml vial Scotchbond™ Universal Adhesive; Instructions for Use; 10 – Sof-Lex™ Spiral Finishing Wheels; 10 – Sof-Lex™ Spiral Polishing Wheels; 1 – RA Mandrel; Technique Guide; Instructions for Use

Refill	Item #	Item #
SHADE	SYRINGE	CAPSULE
A1	038-7126	038-7068
A2	038-7134	038-7076
A3	038-7142	038-7084
B1	038-7159	038-7092
C2	038-7167	038-7100

SYRINGE SHADE REFILLS include 1 – 4g syringe; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS include 20 – 0.2g capsules; Technique Guide; Instructions for Use.

Filtek™ Bulk Fill Flowable Restorative

Refill	Item #	Item #
SHADE	SYRINGE	CAPSULE
A1	502-3767	502-3916
A2	502-3775	502-3924
A3	502-3783	502-3932
Universal	502-3791	502-3940
Item #	Product Description	
502-3890	Dispensing Tips 20 – 19-gauge dispensing tips	

SYRINGE SHADE REFILLS include 2 – 2g syringes; 20 dispensing tips; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS include 15 – 0.2g capsules; Technique Guide; Instructions for Use.

Customer Care Center: 1-800-634-2249 www.3MESPE.com/Filtek

★★★★★ Review these products at www.3MESPE.com/review

3M, ESPE, Filtek, Scotchbond and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2015. All rights reserved.

3M ESPE

Available Through:


PATTERSON
DENTAL

For more information or to place your order, please contact your Patterson representative, local branch or call 1-800-873-7683.