

ESSENTIAL WRIST BLOOD PRESSURE MONITOR



SPEC SHEET

General Description

A classic wrist blood pressure monitor with the benefits of speed and accuracy at the touch of a button.

Product Codes

BV Medical
UB-525

Packaging

Unit Packaging
Single unit, instruction manual, quick start card, and 2 x AAA batteries (included).

Unit Packaging
20 units.

Specifications

Measurement Method: Oscillometric
Measurement Range: Pressure (0-2990 mmHg), Pulse Rate (40-180 beats/min)
Accuracy: Pressure (+/- 3 mmHg), Pulse Rate (+/- 5%)
Power Supply: 2 x AAA batteries (included)
Battery Life: approximately 250 measurements
Warranty: 5 Years
Product Includes: Monitor, Instruction Manual, Quick Start Card, and 2 x AAA batteries

Features & Benefits

QUICK AND EASY MEASUREMENT

Easy one button operation for quick and accurate measurements. Stores 60 measurements in memory with recall. Displays average of all measurements.

TRUSTED ACCURACY

Clinically validated for accuracy. Indicates an irregular heartbeat (IHB) during a measurement. Pressure Rating Indicator classifies blood pressure in accordance with American Heart Association (AHA) medically approved standards.

Barcodes



BV MEDICAL

28W206 COMMERCIAL AVE. UNIT B
LAKE BARRINGTON, IL 60010
(888) 822-8293
FAX: (847) 382-7790
INFO@BVMEDICAL.COM
WWW.BVMEDICAL.COM