

# ENDO 6:1

Endodontics without file clamps

- 
- [HIGHLIGHTS](#)
- [PROPERTIES](#)
- [TECHNICAL DATA](#)
- [DOWNLOADS](#)

## Technical data

Weight	70 g
Reduction ratio	6 : 1
Max. working speed	6,600 min <sup>-1*</sup>
Head height without burr	11.4 mm
Head diameter	8.7 mm
Burr instruments	Ø 2.35 mm RA
Clamping system	Pushbutton chuck
Connection type	ISO interface
REF	59 94 277

\* In combination with a micromotor or at motor speed of 40,000 rpm.