

Why should I switch to GC FujiCEM® Evolve from GC FujiCEM® 2?

GC FujiCEM® Evolve builds on the legacy of GC FujiCEM® 2 which has been used for cementation of more than 17 million restorations¹. It's hard to believe that GC can make a better RMGI, but yes, it has been done! GC FujiCEM® Evolve has all the best properties of the previous generations and more!

¹Based on sold units

Ergonomic syringe delivery (does not require a dispenser):



GC FujiCEM® Evolve



GC FujiCEM® 2

Property comparison:

	GC FujiCEM® Evolve	GC FujiCEM® 2	Difference
Tack cure	Yes	No	Easier clean-up
Shade	Universal	Light yellow	Brighter
Bond strength to Zirconia after thermocycling (MPa)	13.4 (4.5)	6.4 (3.9)	109% higher
Bond strength to lithium disilicate (MPa)	14.0 (4.2)	8.0 (2.9)	75% higher
Radiopacity (AI%)	258	118	119% higher
Flexural strength (MPa)	49 (1)	34 (2)	44% higher

Source: GCC R&D

SKU Comparison for reference:

SKU #	Description	SKU #	Description
012948	GC FujiCEM® Evolve Single Pack	441001	GC FujiCEM® 2 Starter Kit
012949	GC FujiCEM® Evolve Triple Pack Handmix	441002	GC FujiCEM® 2 Refill (Automix Compatible)
012950	GC FujiCEM® Evolve Triple Pack Automix	441000	GC FujiCEM® 2 Automix (42 mixing tips)
012954	GC Push and Click Tip Regular	431578	GC FujiCEM 2 Mixing Tips SL 20/bag
012955	GC Push and Click Tip Endo	431577	GC FujiCEM 2 Mixing Tips Intraoral SL 10/bag