

NEXZR<sup>®</sup>  
TOOL

sagemax



**PERFECT.  
MILLING.**

High edge sharpness  
more details

Long tool lifetime  
Cost reduction and comfort



# NexxZr Tool perfect. milling.

## NexxZr Tool U (Uncoated)

- Solid carbide milling tool
- Double-edge radius cutter
- Optimized cutting areas
- Suitable for ZrO<sub>2</sub>, Wax, PMMA

ZR02

WAX

PMMA



## NexxZr Tool D (Diamond Coated)

- Carbide cutter with special coating for processing of highly abrasive ZrO<sub>2</sub> materials
- Layer thickness approx. 5µm (not electrically conductive)
- Double-edge radius cutter
- Very long tool lifetime

ZR02



NexxZr Tool is a powerful portfolio of dental milling tools. The tool geometries are perfectly matched to the NexxZr zirconia portfolio and enable longer tool lifetime and precise milling results.

The portfolio consists of NexxZr Tool U with an uncoated surface and NexxZr Tool D with a special diamond coating for extremely long tool lifetime.

### NexxZr Tool U

Tool lifetime  
Edge sharpness



### NexxZr Tool D

Tool lifetime  
Edge sharpness



**High edge sharpness**  
**more details**

**Long tool lifetime**  
**Cost reduction and comfort**

for AmannGirrbach



for imes-icore

for Roland

for vhf

for Zenotec

for Zirkonzahn



**Sales & Technical Support EMEA**  
EC Representative AB Ardent

Generatorgatan 8  
195 60 Arlandastad  
Sweden

✉ +46 8 594 412 57  
📧 info@sagemax.com

**sagemax.com**